


力明33歲了, 深受大眾喜愛的 **Gilson Neo** 省力型微量吸管
 即日起至 3/31 **八爪驚喜價 18,800 原價 25,000**

12爪驚喜價 22,800 原價 29,000

(本促銷活動數量有限, 良機稍縱即逝, 敬請把握)

請問哪一家的Pipetman平均可使用 **20年以上** ???



省力再升級, 最多可達50%以上



下半端可高溫高壓滅菌



所有Pipetman皆由五位數以上天平做校正



GILSON M 電動微量吸管

Gilson Test Report

IDENTIFICATION		VOLUMETRIC DATA	
Model:	PIPETMAN M P8X200M	Gasometric results are converted to volumetric data using the corresponding correction factor Z and taking into account the evaporation if applicable:	
Usable volume range:	20 µL to 200 µL	Volume(µL)	10% 50% Nominal
Serial No.:	GL14604	Mean:	20.022 µL 99.917 µL 199.961 µL
Tips:	Diamond D200	Systematic error:	0.022 µL -0.083 µL -0.139 µL
TEST CONDITIONS		Random error:	0.039 µL 0.115 µL 0.123 µL
Balance mF:	39	TEST RESULTS	
Sensitivity:	10 µg	10% volume:	20 µL - PASS
Correction factor Z (µL/mg):	1.0032	50% volume:	100 µL - PASS
Correction for evaporation:	0.035 µL	Nominal volume:	200 µL - PASS
Basis of adjustment:	E2	The PASS status indicates that this pipette performs within Gilson's specifications and consequently its performance is within the maximum permissible errors allowed by ISO 8655-2. Please refer to the user's guide for details.	
Air temperature:	21.0°C		
Barometric pressure:	101.5kPa		
Relative humidity:	60 %		
CONTROL			
Date:	4 OCT 2011		
Inspected by:	BO		


Pipetman出廠時皆通過嚴格的校正測試, 隨貨附出廠報告, 以確保其準確度



GILSON Concept 電動微量吸管



所有Gilson零件、耗材等會印上Gilson Logo及批號, 以確保產品為法國原廠製作, 並符合可追溯的特性



GILSON Trackman 樣品配置導引器



您所熟悉的**力明**, 您所信賴的 **Gilson** (02)2934-2812 Web: www.lmico.com.tw

新竹正明: (03)540-8448 南部利明: (06)267-4330 外縣市專線: 0800-001-088

