

# 最特級的 *Pipetman L*


- L系列維持Pipetman®一貫的高精準度、高品質, 為市面上吸排力道最省力之微量吸管
- 剔除器方向可調整, 並減少其施力
- 專利體積鎖定系統, 確保使用者舒適與可靠的液劑處理



## Patented volume locking system



## Lowest pipetting force on the market



## Trusted results

Pipette model	Volume tested	Maximum Permissible Errors			
		Gilson		ISO 8655	
		Systematic error (μL)	Random error (μL)	Systematic error (μL)	Random error (μL)
P2L	Min.: 0.2	±0.024	≤ 0.012	±0.08	≤ 0.04
	0.5	±0.025	≤ 0.012	±0.08	≤ 0.04
	1	±0.027	≤ 0.013	±0.08	≤ 0.04
	Max.: 2	±0.030	≤ 0.014	±0.08	≤ 0.04
P10L	Min.: 1	±0.025	≤ 0.012	±0.12	≤ 0.08
	5	±0.075	≤ 0.030	±0.12	≤ 0.08
	Max.: 10	±0.100	≤ 0.040	±0.12	≤ 0.08
P20L	Min.: 2	±0.10	≤ 0.03	±0.20	≤ 0.10
	10	±0.10	≤ 0.05	±0.20	≤ 0.10
	Max.: 20	±0.20	≤ 0.06	±0.20	≤ 0.10
P100L	Min.: 10	±0.35	≤ 0.10	±0.80	≤ 0.30
	50	±0.40	≤ 0.12	±0.80	≤ 0.30
	Max.: 100	±0.80	≤ 0.15	±0.80	≤ 0.30
P200L	Min.: 20	±0.50	≤ 0.20	±1.60	≤ 0.60
	100	±0.80	≤ 0.25	±1.60	≤ 0.60
	Max.: 200	±1.60	≤ 0.30	±1.60	≤ 0.60
P1000L	Min.: 100	±3	≤ 0.6	±8	≤ 3.0
	500	±4	≤ 1.0	±8	≤ 3.0
	Max.: 1000	±8	≤ 1.5	±8	≤ 3.0
P5000L	Min.: 500	±12	≤ 3	±40	≤ 15
	2500	±15	≤ 5	±40	≤ 15
	Max.: 5000	±30	≤ 8	±40	≤ 15
P10mL	Min.: 1000	±30	≤ 6	±60	≤ 30
	5000	±40	≤ 10	±60	≤ 30
	Max.: 10000	±60	≤ 16	±60	≤ 30



您所熟悉的**力明**, 您所信賴的 *Gilson* (02)2934-2812 Web: [www.lmico.com.tw](http://www.lmico.com.tw)

新竹正明: (03)540-8448 南部利明: (06)267-4330 外縣市專線: 0800-001-088

